




GUTH LABORATORIES, INC.

590 NORTH 67th STREET • HARRISBURG, PA 17111- 4511 • TELEPHONE: 717-564-5470

Dry Gas Ethanol Breath Standards



Guth Laboratories, Inc. now offers the industry an exceptional dry gas standard. With Guth Laboratories 30 plus years of experience in the breath testing field and the expertise of our gas suppliers in the specialty gas market we are able to offer an exceptional product manufactured in an ISO accredited gas production facility. All standards are NIST traceable and approved by NHTSA.

Gas mixtures are available to support all alcohol breath testers. Contact us for assistance in choosing the correct product for your application.

- NIST National Institute of Standards and Technology traceable
- NHTSA approved
- Available in the following stocked concentrations .038, .040, .080, .082, .100, Inquire for custom concentrations.
- 2 year expiration date
- +/- 2% accuracy
- Full line of regulators and recycling tools available
- Regulators sold separately



Cylinder size	Height	Diameter	Weight	Number of tests	Expiration date from date of manufacture	List Price	Industry standard valve used	Refillable
105 Liter	13.75 inches	3.25 inches	3.5 lbs	300	3 years	inquire	Yes	No
34 Liter	10.25 inches	2.75 inches	1 lb	100	3 years	inquire	Yes	No

**FOR MORE INFORMATION
CALL TOLL FREE 800-233-2338
OR VISIT OUR WEB SITE WWW.GUTHLABS.COM**